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1 Introduction

Please read these instructions carefully. This document is subject to continuous revision. Please verify that this printed version is identical to the current one at ifu.fisso.com. Improper handling can cause harm to the patient or damage the products. Use universal precautions when handling contaminated or biohazardous components/materials.

Manufacturer:

Baitella AG
Thurgauerstrasse 70
CH-8050 Zurich
Switzerland

Phone: +41 44 305 80 00
E-Mail: info@baitella.com
Internet: www.fisso.com

Baitella GmbH
Moltkestrasse 41
DE-78532 Tuttlingen
Germany

UKRP:

confinis (UK) RP Ltd
St John's Innovation Centre
Cowley Road
Cambridge CB4 0WS
United Kingdom
E: ukrp@confinis.com
www.confinis.com

Prescription Device – Rx only

Federal (United States) law restricts these devices to sale by or on the order of a physician.

Notice regarding serious incidents:

For a patient/user/third party in the European Union and countries with an identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national authority.

Copyrights:

All rights reserved. Duplication, adaptation, or translation, of any part of this document, without previous written authorization by Baitella AG is prohibited, except within the framework of the copyright regulations.

Reserved rights to technical changes!

The illustrations and specifications in this instruction for use may differ slightly from actual products.

1.1 Symbols used

Symbol	Definition	Danger	Consequence
	DANGER!	Immediate danger for people	Death or serious injuries
	WARNING!	Possible danger for people or objects	Health damages or severe material damages
	ATTENTION!	Possible danger to objects	Material damages

Symbol	Description
	Symbol for products that are put into use according to the regulation (EU) 2017/745 on medical devices
	Symbol for products that are put into use according to the regulation (UK) UK MDR 2002 on medical devices
	Note Additional hints or important information

	See instruction for use
Rx only	Usage only by trained medical personnel
	Manufacturer
	The articulated arm must not be sprayed directly with liquid or immersed in liquids. The central handle of the articulated arm must be tightened during cleaning/disinfection.
EC REP	European authorized representative
	Distributor
	Importer
MD	Medical device (in this instruction for use indicated as «product»)
UDI	Unique Device Identification
REF	Article number
	Date of manufacture
LOT	Lot number
SN	Serial number
QTY	Quantity
	Keep dry
	Keep away from sunlight
	Humidity limitation
	Temperature limit
	Do not use if the package is damaged

2 Intended Use

This product is used to hold and position breathing circuits, pressure transducers plates and other devices. It must only be used by trained medical personnel capable of judging and controlling any danger to patients. If not the case, the user assumes all responsibility.

	DANGER! Unauthorized changes or modifications to the FISSO holding system are prohibited for safety reasons.
	DANGER! All FISSO components are designed for optimal performance when used together with the FISSO holding system. If a product/component from another manufacturer is used, the user takes full responsibility.
	DANGER! The product will be delivered non-sterile and is non-sterilizable. Prior to the initial and each following use, the product has to be cleaned and disinfected as well as checked for visual irregularities and malfunctioning according to the indications given in this manual.
	WARNING! Do not attach height sensitive devices or transducers with drainage capacity. Injury to patient may occur.
	Guideline for medical devices (conformity declaration) This product corresponds to the regulation (EU) 2017/745 on medical devices.
	Guideline for medical devices (conformity declaration) This product corresponds to the regulation (UK) UK MDR 2002 on medical devices
Rx only	US Federal Law restricts the products to sale by or on the order of a physician.

Special care should be taken when these products are used together with high-frequency applications. Contacts between the product and the high-frequency device must be avoided.


3 Product Description

The following sections will describe the products based on examples. The product range along with all related information can be reviewed on www.fisso.com.

3.1 Holding System (Product)


A complete holding system (Anesthesia) consists of:

- head component (Pos. 1)
- articulated arm (Pos. 2)
- column (Pos. 3)
- rails clamp (Pos. 4)


A complete holding system (Ventilation) consists of:

- head component (Pos. 1)
- articulated arm (Pos. 2)
- column (Pos. 3)
- rails clamp (Pos. 4)

A complete holding system (Monitoring) consists of:

- head component (Pos. 1)
- articulated arm (Pos. 2)
- column (Pos. 3)
- rails clamp (Pos. 4)

4 Installation and use

4.1 Receiving Inspection

Check the product immediately after receipt for eventual transport damages and completeness. Complaints can only be considered if the seller or freight forwarder is immediately notified. In this case, there must immediately send a damage protocol to the next FISSO representative or the company Baitella AG.

WARNING!

Do not use when packaging is damaged.

4.2 Assembly


WARNING!

All handling must be done with as little force as possible and as much force as necessary!

WARNING!

The assembly and handling of the products must be done manually and without additional tools.

1. Attach the rails clamp to the side rail of the surgical table or patient bed and secure it by turning the handle (Pos. 1) from below clockwise. Pay attention during attachment that the upper hook (Pos. 2) grabs behind the rail. Check if the rails clamp is fixed securely and firmly on the side rail.

Step 2 and 3 refer to rails clamps with height-adjustable columns

2. If necessary remove the stop lid (Pos. 3) from the column (Pos. 4) by turning counterclockwise. Place the column (Pos. 4) with one hand into the column opening of the rail clamp and secure the column (Pos. 4) with the other hand by turning the handle (Pos. 2) clockwise. Then screw the stop lid (Pos. 3) into the column (Pos. 4) clockwise and tighten. Check if the column (Pos. 4) is fixed securely and firmly into the rail clamp.
3. The desired height and orientation of the column (Pos. 4) including the articulated arm can be adjusted at any time by loosening the handle (Pos. 5). **Important:** Hold the column (Pos. 4) with one hand and loosen the handle (Pos. 5) with the other hand. When the desired position is reached, tighten the handle (Pos. 5).

WARNING!

If the rails clamp or the column is not correctly fixed, these components may become loose and may cause injuries to the patient.

ATTENTION!

The handle (Pos. 5) should not be tightened at any time without the column first being in place.

4.3 Disassembly

The disassembly of the Arm support must be performed also without additional tools and occurs in the reverse order of the assembly:

Hold the holding system with one hand and loosen with the other hand the fastening screw of the rails clamp. Remove the rails clamp from the rail of the surgical table or patient bed.

4.4 Handling of the Head Components


The device (e.g. breathing circuit, pressure transducer) to be held by the head component must be secured safely.

DANGER!


If the device is not fixed properly it can come loose and move, which can cause injuries. Check secure fixation of the device.

Handling of swinging head

By opening and closing of the wing-screw (Pos. 1) of the swinging head, the movability of the tube holder can be increased or decreased.

4.5 Handling of the Articulated Arms

Holding-Releasing-Positioning-Fixing-Controlling:

1. Hold the articulated arm with one hand in the anterior segment (1) and manipulate the central handle (Z) with the other hand.
2. To release turn the central handle (Z) counterclockwise (2) as far as necessary.
3. Move the articulated arm into the desired position.
4. To secure, turn the central handle (Z) clockwise (3).
5. Controlling: Check that the articulated arm is firmly tightened and functioning properly.

DANGER!

The articulated arm can move unintentionally if the central clamping is released. Always hold the head component on the anterior segment (1) with one hand and manipulate the central clamping device (Z) with the other hand. Always hold the arm support with one hand and manipulate the central clamping device (Z) with the other hand.

DANGER!

If the articulated arm is not tightened properly it can become loose and move, potentially causing injuries. The articulated arm must be tightened properly.

DANGER!

A damaged product can lead to serious injuries. Only use products in perfect condition and check functionality.

WARNING!

Keep away your extremities from the joints during tightening, potentially causing injuries. Keep away your extremities from the joints during untightening and transportation, potentially causing injuries.

WARNING!

The articulated arm can transmit electrical current and heat to the patient. Avoid the contact between the articulated arm and any sources of electrical current or heat.

ATTENTION!

The fixation of the articulated arm is based on the principle of friction. Changing the position without loosening the clamping mechanism can cause damage and will shorten the lifespan of the articulated arm. The articulated arm can be adjusted with little force. If the central clamping handle (Z) has been fully loosened it has to be turned in the clockwise direction!

ATTENTION!

If the articulated arm is fixed to an apparatus, the apparatus must not be moved/shifted by means of the articulated arm!

Do not use the articulated arm to move/shift an apparatus.

5 Processing

All devices are to be cleaned and disinfected before each application; this is required as well for the first use after delivery of the device (cleaning and disinfecting after removal of the protective packaging). Effective cleaning and disinfection is an indispensable requirement for the safe application of the products.

You are responsible for the safe application of the device. Therefore, please ensure that only sufficiently product-specifically validated procedures will be used for cleaning and disinfection as well as that the validated parameters will be applied for each cycle.

Additionally, please pay attention to the legal provisions valid for your country as well as to the hygienic instructions of the hospital or doctor's practice.

	<p>DANGER! The product will be delivered non-sterile and is non-sterilizable. Before the initial and each following use, the product has to be cleaned and disinfected as well as checked for visual irregularities and malfunctioning according to the indications given in this manual.</p>
	<p>ATTENTION! Material damage due to improper cleaning and disinfection. This product can be cleaned and disinfected only manually. The cleaning and disinfection agents must be compatible with the product.</p> <p>Do not expose any product to temperatures higher than 55°C (131 °F).</p>
	<p>WARNING! Danger of explosion. Alcohol-containing agents form ignitable mixtures, which can explode if used together with high-frequency applications. Do not use alcohol-containing agents together with high-frequency applications.</p>
	<p>WARNING! The articulated arm must not be sprayed directly with liquid or immersed in liquids. The central handle of the articulated arm must be tightened during cleaning/disinfection.</p>

Material resistance

When choosing cleaning and disinfecting agents, please make sure that they do not contain the following:

- organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- strong lyes (maximum admitted pH-value 12, neutral/enzymatic cleaning detergent recommended)
- organic solvents (for example: acetone, ether, alcohol, benzine)
- oxidizing agents (for example: peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic, halogenated hydrocarbons
- oils

Please consider during the selection of the detergents in addition, that any ingredients deviating from ethanol and isopropanol may cause potential critical remnants on the device.

Check

Check after cleaning, respectively cleaning/disinfection the functionality of all components of the product to ensure smooth operation throughout the intended range of motion, corrosion, damaged surfaces, chipping, and dirt, and sort out damaged products (limitation of re-use, see section 5.3 Re-usability/ End of Life). Products that are still dirty must be cleaned and disinfected again.

Sterilization

Any sterilization procedures must not be applied.

Packaging

Please cover the cleaned and disinfected device with single-use protective bags, which fulfill the following requirements:

- foil bag

- made of plasticizer-poor plastic material and with suitable cleanliness (at least food grade)

Label the bag with the cleaning and disinfection status.

Packaging is not required if the device is immediately used.

5.1 Cleaning and disinfection

5.1.1 Pre-treatment

Please remove coarse impurities of the products directly after use (within a maximum of 2 hrs). To this end, wipe all visually contaminated surfaces of the products intensively with disinfection wipes¹.

¹ Please consider that the wipes should be alcohol- and aldehyde-free (otherwise there is a risk of fixation of blood impurities or damage to the product), possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), be suitable for disinfecting instruments made of metallic or plastic material and compatible with the products (see chapter "Material resistance"). Please consider that a disinfectant used in the pre-treatment step is solely for the safety of the staff and cannot replace the disinfection step later to be performed after cleaning.

5.1.2 Manual cleaning and disinfection by wiping

Observe the following points when selecting the wipes for cleaning and disinfection:

Fundamental suitability for the cleaning and disinfection of products made of metallic or plastic material
Approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking)
Compatible with the products (see chapter Material resistance)

NOTE

The instructions for use of the cleaning and disinfecting agent manufacturer must be followed (including instructions regarding concentration, temperature, and soaking time as well as post-rinsing).

Proof for the general suitability of the device for an effective cleaning and disinfection was provided by an independent test laboratory using the cleaning and disinfecting mikrozid sensitive wipes (Schülke & Mayr GmbH, Norderstedt) was applied.

Procedure:

1. The central handle of the articulated arm must be tightened during cleaning/disinfection
2. Wipe all surfaces of the product completely and carefully with fresh wipes. Visibly dirty or dry wipes must not be used further.
3. Check all surfaces for remaining visible contamination and wipe again by applying a fresh wipe, if required.
4. Check the products (see chapters "Control" and "Maintenance").
5. Reassemble the product and wipe again using a fresh wipe.

5.1.3 Manual cleaning and disinfection – for tube holder (ventilation)

Observe the following points when selecting the cleaning and disinfection detergents:

- Fundamental suitability for the cleaning and disinfection of instruments made of plastic material
- Disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible with the cleaning detergent used
- Compatibility of the used detergents with the products (see chapter "Material resistance")

Combined cleaning/disinfection detergents are not recommended.

NOTE

The instructions for use of the cleaning and disinfecting agent manufacturer must be followed (including instructions regarding concentration, temperature, and soaking time as well as post-rinsing)..

Procedure:

1. Cleaning
 - a. Remove the tube holder from the holding system.
 - b. Soak the tube holder for the given soaking time in the cleaning solution so that the tube holder is sufficiently covered. Assist by careful brushing with a soft brush (never use metal brushes or steel wool), ultrasonic treatment is not recommended. Ensure that there is no contact with other products.
 - c. Then remove the tube holder from the cleaning solution and post-rinse it at least three times with water.
 - d. Check the tube holder (see chapters “Control” and “Maintenance”).
2. Disinfection
 - a. Soak the tube holder for the given soaking time in the disinfectant solution so that the tube holder is sufficiently covered. Ensure that there is no contact with other products.
 - b. Then remove the tube holder from the disinfectant solution and post-rinse it at least five times with water.
 - c. Dry and pack the tube holder immediately after removal (see chapter “Packaging”).

5.2 Storage

WARNING!

Store the products such that they are not damaged.

Please store the device after cleaning, disinfection, and complete drying in the protective packaging at a dry and dust-free place.

These storage conditions are not required if the device is stored directly at the location of the application and will be used immediately.

5.3 Re-usability / End of life

The products can be reused for 10 years if they are not damaged and have been cleaned and disinfected according to the instructions provided. Every re-use beyond 10 years and the use of a damaged or dirty product are the responsibility of the user. If a damaged product is used or if a product is re-used without having been cleaned and disinfected, then the company is not liable for any damage caused.

In case of disrespect, all liability is excluded.

6 Storage, Handling, and Transportation

WARNING!

Store the products such that they are not damaged.

The product must be stored and transported in dry conditions in the original packaging at room temperature and not exposed to direct sunlight. Incorrect storage and transportation may influence the product characteristics leading to failure.

7 Maintenance

Special maintenance is not required. The product should be sent back to the manufacturer or distributor if the action of moving parts or the clamping force is impaired or in case of damage.

Instrument oil or grease may never be used.

NOTE

Repair must be carried out by Baitella AG. If not respected, the warranty becomes void!

ATTENTION!

Please send back only processed products to the distributor/manufacturer for repair (cleaned and disinfected).

8 Disposal

The products must be correctly disposed of under national regulations and medical guidelines.

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All rights reserved. Baitella AG 의 사전 서면 승인 없이는, 명확하게 저작권 규제에 저촉되지 않는 경우를 제외하고 이 문서의 어떠한 부분도 복제, 각색 또는 변환/번역해서는 안 됩니다.